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Literatur

1. Magnusson, S., Baldursson, B. T., Kjartansson, H., Rolfsson, O. & Sigurjonsson, G. F. Regenerative and Antibacterial Properties of Acellular Fish Skin Grafts and Human Amnion/Chorion Membrane: Implications for Tissue Preservation in Combat Casualty Care. *Mil. Med.* 182, 383–388 (2017).
2. Magnusson, S. et al. Decellularized fish skin: characteristics that support tissue repair. *Laeknabladid* 101, 567–573 (2015).
3. Rakers, S. et al. Fish matters: the relevance of fish skin biology to investigative dermatology. *Exp. Dermatol.* 19, 313–324 (2010).
4. Baldursson, B. T. et al. Healing rate and autoimmune safety of full-thickness wounds treated with fish skin acellular dermal matrix versus porcine small-intestine submucosa: a noninferiority study. *Int. J. Low. Extrem. Wounds* 14, (2015).
5. Kirsner, R. S. et al. Double-Blind, Prospective, Randomized Clinical Trial on 170 Acute Wounds Shows Significantly Faster Healing Rate with Intact Fish Skin Compared to Human Amniotic Membrane. *Natl. Am. Podiatr. Med. Assoc. Annu. Sci. Meet.* (2018).
6. Lullove E. J. et al. A multicenter, blinded, randomized controlled clinical trial evaluating the effect of Omega-3-rich fish skin in the treatment of chronic, nonresponsive diabetic foot ulcers. *Wounds*. Published online April 15, 2021.
7. Winters C, Kirsner RS, Margolis DJ, Lantis JC. Cost Effectiveness of Fish Skin Grafts Versus Standard of Care on Wound Healing of Chronic Diabetic Foot Ulcers: A Retrospective Comparative Cohort Study. *Wounds*. 2020;32(10):283-290.
8. Stone R 2nd, Saathoff EC, Larson DA, et al. Accelerated Wound Closure of Deep Partial Thickness Burns with Acellular Fish Skin Graft. *Int J Mol Sci.* 2021;22(4):1590.
9. Pujji O, Jeffery SLA. Safe burn excision prior to military repatriation: an achievable goal? *BMJ Military Health* 2018;164:358-359.
10. Chun K, Yang, John C, Lantis II & Thais O, Polanco. A prospective, single-center, non-blinded, non-comparative, post-market compassionate clinical evaluation of a Novel Acellular Fish Skin Graft which contains Omega3 fatty acids, for the closure of hard to heal lower extremity chronic ulcers. *Wounds* 28, 112–118 (2016).
11. T. T. Trinh, F. Dünschede, C.-F. Vahl & B. Dorweiler. Marine Omega3 Wound Matrix for the Treatment of Complicated Wounds. *Phlebologie* 45, 93–98 (2016).
12. Dorweiler, B. et al. Die marine Omega-3-Wundmatrix zur Behandlung komplizierter Wunden. *Gefässchirurgie* 22, 558–567 (2017).
13. Woodrow, T., Chant, T. & Chant, H. Treatment of diabetic foot wounds with acellular fish skin graft rich in omega-3: a prospective evaluation. *J. Wound Care* 28, 76–80 (2019).
14. Sitje, T. S., Grøndahl, E. C. & Sørensen, J. A. Clinical innovation: fish-derived wound product for cutaneous wounds. *Wounds Int.* 2018 9, 44–50 (2018).
15. Patel, M. & Lantis II, J. C. Fish skin acellular dermal matrix: potential in the treatment of chronic wounds. *Chronic Wound Care Manag. Res.* 6, 59–70 (2019).
16. Sibbald R, Goodman L, Woo K. et al (2011) Special considerations in wound bed preparation 2011: An update. *Adv Skin Wound Care.* 24(9):415-36.
17. Wallner C. et al. A Comparison of Intact Piscine Skin, Split-thickness Skin Graft, and Lactic Acid Membrane in Treating Superficial and Deep Burn Wounds Following Enzymatic Debridement. *J Burn Care Res.* 2021; 42 (Suppl 1): 125- 126

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