

EMPOWER
life changing decisions.



DxFLEX – System Overview


Beckman Coulter Life Sciences, Flow Cytometry

DxFLEX flow Cytometer is CE marked for up to 13-color in vitro diagnostic use. This device is not available in all countries. Please check with your local sales representatives before placing your orders.

FLOW-7996CP10.20

Empowering those seeking answers to life's important scientific and healthcare questions.

Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks used herein are the trademarks or registered trademarks of Beckman Coulter, Inc. in the United States & other countries.



DxFLEX E-Roadshow, Oct 2020





DxFLEX: Additional Features (over CytoFLEX*)

- Key features:
 - CE-marked for clinical use for up to 13 “colors”
 - 4 upgradable configurations: 1L5C, 2L6C, 2L9C, 3L13C
 - Optional Autoloader for 32 position MCL carousel
- New clinical **CytExpert for DxFLEX software**:
 - User and account management
 - Experiment and report templates
 - Panel experiments and flagging of out-of-range results
 - ...



Compensation Library Manager

Fluorescence: FITC

Compensation Keywords: TEST 2014-09-18 12:38:32

Label: Lot No:

Keywords: TEST Gain: 293

Cell Autofluorescence

Height: 0.00 Area: 0.00

Bead Autofluorescence

Height: 1.60 Area: 1.02

Fluorescence	Height	Area	Gain
PE	11.15 %	11.68 %	152
PerCP-Cy5.5	2.09 %	2.22 %	687
PE-Cy7	0.14 %	0.22 %	742
APC	0.12 %	0.21 %	352
APC-H7	0.11 %	0.21 %	758
APC-cf7	0.01 %	0.03 %	500
BV421	0.00 %	0.00 %	289
BV510	1.44 %	1.95 %	275
BV570	0.11 %	0.24 %	222
BV570	0.07 %	0.06 %	1287
BV711	0.02 %	0.01 %	921

Export... Import... Delete OK Cancel Apply

DxFLEX is based on the highly successful CytoFLEX* technology, but includes additional features needed for clinical work



Empowering those seeking answers to life's important scientific and healthcare questions.

Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service name used herein are the trademarks or registered trademarks of Beckman Coulter, Inc. in the United States & other countries.

DxFLEX is not available in all countries. Not available in the US.

* For Research Use Only. Not for use in diagnostic procedures.






What is different about DxFLEX compared to other clinical flow devices available on the market?

- **Wavelength Division Multiplexer (WDM)**


- Compact design with short travel distance of fluorescence detection between the detectors
- High yield light transmission concept without dichroic filters



Bandpass only for flexibility and efficiency




Mirror band pass filters
instead of
dichroic/bandpass pairs



FLOW-4067SB09.18



Ad 1: Better sensitivity in the far red fluorescence channels



Why are APDs better than PMTs?

Because the quantum efficiency of APDs in the spectrum of commercial available dyes is higher

This is especially true for higher wavelength.

Image comes from "Introduction to Beckman Coulter's CytoFLEX* Research Flow Cytometer Platform" presentation. FLOW-3070CP09.17
* For Research Use Only. Not for use in diagnostic procedures





Empowering those seeking answers to life's important scientific and healthcare questions.

Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service mark used herein are the trademarks or registered trademarks of Beckman Coulter, Inc.



Ad 1: What does it mean for your daily work?

Better sensitivity to analyze dim populations in all fluorescence channels



■ Blue – APD (DxFLEX)
■ Red – PMT (Navios)

Rainbow 8 peak beads analysis: FL5 (Blue Laser)


Empowering those seeking answers to life's important scientific and healthcare questions.

Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks used herein are the trademarks or registered trademarks of Beckman Coulter, Inc. in the United States & other countries.



Ad 1: What does it mean for your daily work?

Better sensitivity to analyze dim populations in all fluorescence channels





Rainbow 8 peak beads analysis: FL8 (Red Laser)

Empowering those seeking answers to life's important scientific and healthcare questions.

Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service mark used herein are the trademarks or registered trademarks of Beckman Coulter, Inc. in the United States & other countries.



Ad 1: Sensitivity APD vs PMT Technology



Whole blood analysis of the same prepared sample to analyse NK cells

CD56-PC7 vs CD3-PB

■ Blue - DxFLEX
■ Purple - FACSCanto

Better APD sensitivity especially in the far red channels.

Data kindly provided by Dottoressa Bertaina, Hospital Bambin Gesù - Rome

Empowering those seeking answers to life's important scientific and healthcare questions.


Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service mark used herein are the trademarks or registered trademarks of Beckman Coulter, Inc. in the United States & other countries.

* The scales of representation of the axes cannot be set at the same values because the different number of logarithmic scales available for displaying the data by the systems (DxFLEX 7 decades VS FACSCanto 5 decades)





Ad 2: Spreading comparison - DxFLEX Vs PMT device



Data kindly provided by Ospedale Gemelli, Rome

Empowering those seeking answers to life's important scientific and healthcare questions.

Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marked herein are the trademarks or registered trademarks of Beckman Coulter, Inc. in the United States & other countries.


*The spreading effect can be evaluated only from a visual point of view. Because the different number of logarithmic scales available for displaying the data doesn't allow a comparable statistical analysis.








Conclusions


1. Better sensitivity due to the improved APD Quantum Efficiency.
 - Better separation of dim populations especially in the Far Red channels.
 - Less complexity in panel design.
2. Less fluorescence spreading
 - Better resolution of dim fluorescent populations.
3. Linearity response to gain variation
 - Gains can be changed at any time because compensations are automatically updated.
 - Simplified workflow process for compensation setup.



Avalanche Photodiode.
The DxFLEx flow cytometer uses Avalanche Photodiode detectors instead of PMTs. The low electronic noise contributes to the resolution capabilities of the instrument.



Video 1



Video 2